

## Using the Eye Gaze

The Eye Gaze works like a mouse – but using your eyes!

The Eye Gaze Computer is kept in the Den, on height adjustable stand.

Keyboard, mouse and Eye Gaze bar kept in yellow locker in the Den – key in Team Room. Orange label on key.

Plug in Computer, turn on – switch on side of computer.

Check mouse has been turned on.

When 'Rainbow' picture appears, swipe up from bottom of screen to get login screen.

Enter Lifelites password: Life0000

Connect Eye Gaze bar to magnetic strip at bottom of computer, connect lead to lead hanging from side of computer.

Lower computer stand so that Eye Gaze bar is roughly in line with person's head. You can bring someone in nearer to the screen too.

Click on Eye Gaze Icon,



then click on 'Track Status'

A Black box will appear, ask the person using EyeGaze just to keep looking at the screen. Two white dots should appear in the black box, they need to be in the 'green' zone as shown on the side of the black box. Adjust the height of the computer, bring the person nearer to get the two white dots into the right position.

When you have got it adjusted correctly, minimise this screen.

Now click on the different programmes you can use with EyeGaze.

To Start with try:

**Eagle Paint (EP icon)**



Toggle space bar to start the programme. As the person looks at the screen, stars/circles/squares appear, eye movement will make the shapes move across the page.

To change the shapes – go to Drawing Mode and click on shape you require. Click on Background to change from white screen to black or v.v. Click on Edit then clear screen to start again.

To Print: go to File, click on Print, write name in box, click on Preview, then click on print. Should print out of Lifelite printer.



**Jackson Pollock Paint.**

(nb if you use the mouse it will make a 'mark' on the screen.)

Just gaze at the screen and see what happens.

To print: click on print icon.



**Eye Gaze Aliens**

By gazing at the Aliens, they should 'splatt'!

See the booklet How to use EyeGaze with HelpKidzLearn for next level of games, developing concentration and to alter the set up.

Important:

Use Eye Gaze for short bursts of time – 10 – 15mins max as it is very tiring on the eyes.

Give the computer and Eye Gaze system time to get going.

Handy Hint: Hold a piece of paper over the Eye Gaze to stop it working when you need to have a break!

